


GUIDE WIRE




Guide wires are used to gain ureteral access to establish a tract and to assist in the placement, replacement and exchange of medical devices during urological procedures. These guide wires are not intended for PTCA (Percutaneous transluminal coronary angioplasty) use.

Guide wires are available in Stainless Steel PTFE coated Ever Green, Nitinol Hydrophobic coated, Zebra and Black & White.

Guide wires are delivered in sterile peel-open packages. Intended for One-time Use.

GUIDE WIRE



Straight Tip Guide Wire



J Tip Guide Wire

Ordering Details :

Family	Size	Length (cm)	Color / Material
GW	0.018" 0.025" 0.028" 0.032" 0.035" 0.038"	(80 PTFE only) 150 (250 PTFE only)	*N *BW *Z *PEG

J J-Tip
S Straight Tip

*PEG Polyurethane Ever Green
*N Nitinol (Hydrophilic coated)
*BW Black & White
*Z Zebra

Medpro MEDICAL B.V.
Landjuweel 11, 3905PE Veenendaal
the Netherlands
T : +31 318 769080 F : +31 318 769107
E : sales@medpromedical.nl
w : www.medpromedical.nl